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GENERAL DESCRIPTION

The Dial medical external fixation system consists of external fixators-unilateral, circular, hexapod- and PinFix®/PinWire bone trays that acts as a link between the fractured bone and the external fixator allowing forces normally supported by the bone (vertical, lateral, rotational forces) to be transferred to the external fixator to which they are connected by specially designed components secured by nuts and bolts.

External fixation systems are made to stabilize bone fractures and/or perform bone transport and/or angular corrections following post- traumatic or genetic deformities of long bones. The purpose is to mechanically stabilize the bone so that it can achieve healing or be modified in axial or angular elongation consequent to orthopedic treatment according to the technique of callogenesis, provide mechanical support to the bone to stabilize fractures and perform corrections of deformities.

THE EXTERNAL FIXATION SYSTEM BY DIAL MEDICALI CAN BE USED BY ONLY SPECIALIZED MEDICAL DOCTOR.

Transosseous paths PinFix® are available in different configurations, conical and cylindrical, with different diameters (with and without thread), in stainless steel or titanium; the thread can be coated with hydroxyapatite.

Transosseous paths PinFix® are available both self-tapping and self-drilling

For a safe and correct use please refer to the pre-cautions and tips from clinical literature and surgical techniques provided by Dial medicali.

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- Presence of infection
- Patients with mental instability, such that they cannot follow the surgeon's instructions for the postoperative course
- HIV-positive patients
- Patients with severe osteoporosis
- Patients presenting with diabetes mellitus
- Patients who have shown a particular sensitivity or allergy to the PinFix® implant material

IMPORTANT NOTES

- PinFix® has two configurations: cylindrical and conical

- When you use cylindrical pin you can retract it from the bone without losing its seal

- When you use conical pin never retract it after its insertion to avoid mechanical loosening

Otherwise the conical pin can be screwed deeper to increase grip in the bone during treatment.

- Not all surgeries are successful. Further complications can develop at any time due to misuse, medical

reasons, or device failures, resulting in the need for new surgery to replace or remove the device.

Pre-operative and operative procedures, which require knowledge of surgical techniques, the correct choice and positioning of the devices are fundamental factors for the success of the surgery. Correct information to the patient on the post-operative course is an important factor for the success of the operation. Not all surgeries are successful. Further complications can develop at any time due to misuse, medical reasons, or device failure, resulting in the need for reoperation to replace or remove the device. Preoperative, operative, and postoperative procedures that require knowledge of surgical techniques, proper device selection, and device placement are critical factors for successful surgery. Proper patient information about the postoperative course is an important factor in successful surgical treatment.

WARNINGS AND PRECAUTIONS

- In pediatric applications pay attention to the positioning of the screw in the metaepiphyseal area to prevent possible interference with the growth cartilages (physis).
- Select the diameter and length of the screw based on the size of the bone; the PinFix® screw with D. 6 mm should be used on bone with a diameter of at least 18/20 mm, otherwise evaluated to use a smaller caliber screw.
- If the screw is inserted by hand, use dedicated key or spindle-Jacobs type.
- If the screw is inserted by a motor, do not use a high-speed rotation.
- The insertion of the pin-wire must take place by rotation of the same using hand or electric power and not by percussion
- To assure the maximum hold the screw has to always bi-cortex with at least 1-2 threads showing. Not respecting this procedure reduces the mechanical stability of the screw (around 30%-50%). Moreover, it can cause premature osteolysis.
- To prevent osteolysis and to even out the load between the bone and the fixator the screws should be positioned perpendicular to the diaphysis or parallel the articular joint.
- Self-drilling screws do not normally require pre-drilling but the same may be necessary if the cortical bone is particularly thick; if necessary, use 4.5 / 4.8mm perforators for 6mm diameter screws and 3.2 / 3.5mm perforators for 5mm diameter screws; if necessary, screws with a diameter of 3 - 4 mm require a perforator with a diameter of 1.0 / 1.2mm smaller than the gauge of the screw used.
- Stainless steel PinFix® screw contains Nickel. It is recommended that the patient undergoes nickel allergy tests if there is any doubt.
- It is mandatory to inform patients about the hygiene to be carried out at the screw application site.
- When removing the screws, especially if coated with hydroxyapatite, the surgeon is required to evaluate any patient sedation to prevent any discomfort / pain sensation in the patient himself
- Pay particular attention not to penetrate the joint with the screws
- The PinFix®/PinWire medical devices have not been evaluated for safety and compatibility in the magnetic resonance (MR) environment and has not been tested for heating or migration to the MR environment.
- **Following the tests performed, the TFix emergency/urgency fixator has been found to be compatible with the MR environment (including its titanium alloy bone tracts). Other models of external fixators and related intraosseous tracts by Dial Medicali have not been tested for heating or migration in the MR environment.**

RESIDUAL RISKS

- Dislocation, subluxation, joint stiffness
- Fractures caused by the screw holes themselves, following their removal
- Fracture of the screw
- Tissue necrosis
- Postoperative bleeding or vascular disorders
- Damage to nerves or blood vessels, resulting from the insertion of the bone paths
- Bone infection, osteomyelitis or septic arthritis
- Edema or swelling, possible compartment syndrome
- Failure of bone regeneration, development of pseudoarthrosis or nonunion
- Bone deformity
- Fractures of regenerated bone

Complications include the inability to stabilize the extremities. Failure to provide definitive stabilization occurs due to inappropriate pin placement or inadequate structure. It is most commonly seen when there is extensive bone loss and/or loss of ligamentous support around the knee or ankle, in which case additional pins, double stacking, or a multiplanar device may be needed.

The main difficulty that occurs biomechanically is the distance of the fixator from the fracture if the reduction does not provide intrinsic stability. Injury to nerves and vessels can occur due to direct injury to the pin or by traction. Infection of the pin tract correlates with time, and if the pin is left long

enough, there is a high probability of it occurring. With short-term temporary fixation, pin tract infections are less of a concern.

It is suggested that specific texts be consulted for the proper identification of the appropriate anatomic areas for insertion of the intraosseous tramlines.

The Surgeon must always consider the type of fitting he is performing in relation to the treatment for which it is intended, particularly to achieve the best possible compromise of mechanical (and therefore bony) stability with respect to the patient's subjective state and general condition.

The Dial medicali's fixators should be positioned at least 15-20 mm from the skin to prevent overpressures and necrosis.

THE PRODUCTS PROVIDED IN STERILE STATE CANNOT BE RE-STERILIZED.

Users are required to report any serious incidents occurring in relation to the device to the manufacturer and the competent authority of the Member State where the user and / or patient is established.

Materials

The external fixators are made of stainless steel, aluminum alloy, titanium or carbon fiber components.

PinFix® device is made provided in stainless steel 316LVM or Titanium Ti6Al4V (titanium). PinFix®HA has thread with hydroxyapatite (Osprovit® by Lincotek SpA).

PinWire is provided in stainless steel 316LVM or Titanium Ti6Al4V

STERILE AND NOT STERILE PRODUCTS

Devices supplied sterile bear a label indicating this status; unless otherwise specified, devices are NOT STERILE.

Sterile packaging:

External Fixators: Each individual component comes packaged in a double sealed envelope inserted in a box for further protection. Sterilization is by ethylene oxide. Sterility is guaranteed within the stated expiration date and with intact packaging.

PinFix®/PinWire: each individual component comes packaged in a double sealed envelope inserted in a box for further protection. Sterilization is done by irradiation. Sterility is guaranteed within the stated expiration date and with intact packaging.

If in doubt, do not use the product and replace it. **THE PRODUCT CANNOT BE STERILIZED AGAIN AFTER OPENING THE PACKAGE.**

Nonsterile devices: each component comes in dedicated metallic box, cases or plastic bag.

CLEANING, MAINTENANCE AND STERILIZATION

EXTERNAL FIXATORS AND BONE TRACERS IF SUPPLIED NON-STERILE MUST BE STERILIZED BEFORE THEIR USE BY AUTOCLAVING.

Prior to sterilization, products must be carefully examined and washed, then properly packaged for sterilization according to official procedures validated in health care facilities.

CAUTION: Aluminum alloy components are damageable by the use of detergents and alkaline solutions (pH>7). The anodized coating is damageable by cleaning solutions containing free hydrogen ions or sodium hydroxide. It is absolutely not recommended to clean aluminum alloy elements with chloride, iodide, bromide, fluoride-based disinfectants.

Sterilize in steam autoclave at 134C° for 15min.

BOTH THE EXTERNAL FIXATOS AND PinFix®/PinWire ARE EXPRESSLY SINGLE-USE PRODUCTS AND CANNOT BE REUSED. ANY REUSE OF THE DEVICE MAY CAUSE RISK OF INFECTION AND CONTAMINATION. IN ADDITION, THE PHYSICAL AND MECHANICAL CHARACTERISTICS OF THE DEVICES CANNOT BE GUARANTEED.

Dial medicali srl is responsible for safety and effectiveness only when first used.

Follow current regulations for disposal.

Products are considered new when taken from a Dial Medicali original package

Dial medicali's products, depending on how they are provided, they are classified as:

External fixators:

NOT STERILE: Medical device class 1 non-sterile; STERILE: Medical device class 1s sterile.

Transosseous paths (both STERILE and NOT STERILE): medical device class IIb.

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