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INSTRUCTIONS AND PRECAUTIONS FOR USING THE MEDICAL DIAL EXTERNAL FIXATION SYSTEM

MANUFACTURER: DIAL MEDICALI SRL - MILAN

DESCRIPTION AND INTENDED USE

The Dial medicali external fixation system consists of a series of external fixators - unilateral, circular, hexapodal - made to stabilize bone fractures and / or perform bone transport and / or angular corrections following post-traumatic or congenital deformities of the long bones. The Dial medicali external fixation system can only be used by specialized medical personnel.

CONTRAINDICATIONS

The use of Dial medicali external fixation system is contraindicated in patients suffering from mental instability.

WARNINGS AND PRECAUTIONS

Consultation of spec texts is recommended for the correct identification of the appropriate anatomical areas for the insertion of the intraosseous fixations

The surgeon must always consider the type of mounting he is performing according to the treatment he would like to do in particular to achieve the best possible compromise for mechanical (and therefore bone) stability.Dial Medicali Fixators should be placed at least 15-20mm from the skin to prevent overpressure and necrosis.

Dial medical external fixations have not been evaluated for safety and compatibility in a magnetic resonance (MR) environment and have not been tested for heating or migration in the MR environment.

DIAL MEDICALI FIXATORS PROVIDED STERILE MUST NOT BE RE-STERILIZED.

MATERIALS

Dial medicali fixators consist of components in stainless steel, aluminum alloy, titanium or carbon fiber.

STERILE AND NON-STERILE PRODUCTS

The devices supplied sterile have a label indicating this status; unless otherwise specified, the devices are non-sterile.

PACKAGING:

Sterile devices: each single component is supplied packaged in a double bag. Sterilization takes place with ethylene oxide. Sterility is guaranteed within the expiry terms indicated and with the packaging intact. If in doubt, do not use the product and replace it. It cannot be sterilized again after opening the package.

Non-sterile devices: each single component is supplied in dedicated boxes or in a transparent labeled plastic bag.

CLEANING, MAINTENANCE AND STERILIZATION

IF SUPPLIED NON-STERILE THE COMPONENTS OF FIXATORS MUST BE STERILIZED BEFORE USE.

Before sterilization, the products must be carefully examined and washed, then properly packaged for sterilization, according to the official validated procedures. Sterilize in a steam autoclave at 134C ° for 15min.

The aluminum alloy components can be damaged by the use of detergents and alkaline solutions (pH> 7). The anodized coating can be damaged by cleaning solutions containing free hydrogen ions or sodium hydroxide. It is absolutely not recommended to clean the aluminum alloy elements with disinfectants based on chlorides, iodides, bromides, fluorides.

The products are considered new when they are taken from an original Dial Medicali package.

**ALL DIAL MEDICALI FIXATORS ARE DISPOSABLE AND SHOULD NOT BE REUSED.** Dial medicali srl is only responsible for the safety and efficacy in case of first use of the components. For disposal, comply with the effective regulations.

Dial medicali fixators, depending on how they are supplied, are classified as:

Class 1 non-sterile medical device (in compliance with EU Medical Device Regulation 2017/745 on medical devices

Medical device class 1s sterile (in compliance with Directive 93/42 EEC, implemented with Legislative Decree no. 46 of 24/02/97).

**RESIDUAL RISKS AND POSSIBLE COMPLICATIONS**

- Damage to nerves or blood vessels, resulting from the insertion of bone passages
- Bone infection, osteomyelitis or septic arthritis
- Edema or swelling, possible compartment syndrome
- Lack of bone regeneration, development of nonunion or nonunion
- Bone deformity
- Fractures of the regenerated bone

- Inability to stabilize extremities due to inappropriate positioning of the traps or inadequate structure. It is most commonly seen when extensive bone loss and / or loss of ligament support occurs around the knee or ankle.
- Complications include the inability to stabilize the extremities. Failure to provide definitive stabilization occurs due to inappropriate post placement or inadequate frame. It is most commonly seen when extensive bone loss and / or loss of ligament support occurs around the knee or ankle. When this occurs, additional pins, double stacking, or a multiplanar device may be required.

- The main difficulty that occurs biomechanically is the distance of the bar from the fracture if the reduction of the fracture does not offer adequate stability. Injury to the nerves and vessels can occur due to direct pin injury or traction. Pin infection is time-related, and if the pin is left on long enough, there is a high risk it occurs. With short-term temporary fixation, pin infections are rare.

**IMPORTANT:**

Not all surgeries are successful. Further complications can be developed at any time due to incorrect use, for medical reasons or device failures, resulting in the need of a new surgery to replace or remove the external fixation device. Pre-operative and operative procedures, which require knowledge of surgical techniques, the correct choice and positioning of the devices are fundamental factors for the success of the surgery. Correct information to the patient postoperatively is an important factor for the success of the surgery.

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