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PINFIX PATHS FOR EXTERNAL FIXATION

PinFix® is an intrasosseous path also called half-pin, screw, pin-wire or k-wire, which acts as a connection between the fractured bone and the external fixator allowing to transfer the forces normally supported by the bone (vertical, lateral forces , rotational) to the external fixator to which they are connected with specially designed components and fixed by nuts and bolts. The aim is to mechanically stabilize the bone so that it can achieve healing or be modified in axial or angular elongation as a result of an orthopedic treatment according to the callogensis technique, provide mechanical support to the bone to stabilize fractures and perform deformity corrections.

DESCRIPTION AND INSTRUCTIONS FOR USE

The PinFix® devices are a medical device that can only be used with external fixation (circular frame, monolateral and damage control) for the upper and lower limbs, pelvis hand and foot and it can only be used by specialized medical personnel.

PinFix® is available in different configurations, conical and cylindrical, in different diameters (with or without thread) made of stainless steel or titanium and it can have the thread coated with hydroxyapatite.

PinFix® is available in self-drilling and non self-drilling.

For its use refer to the precautions and suggestions reported in the literature and to the informative brochures provided by Dial medicali.

STERILE AND NON-STERILE PRODUCTS

The devices supplied sterile have a label indicating this status; unless otherwise specified, Pinfix devices are non-sterile.

METHOD OF STERILIZATION

PinFix® sterile is packaged in a double sterile bag and an external protective box.

Sterilization takes place by irradiation.

Sterility is guaranteed within the expiry terms indicated and with the packaging intact. If in doubt, do not use the product and replace it. It cannot be sterilized again after opening the package.

Non-sterile PinFix® is supplied either in a box or in a clear plastic bag.

Before sterilization, the products must be carefully examined and washed, then properly packaged for sterilization, according to the official validated procedures. Sterilize in a steam autoclave at 134°C for 15min.

PinFix® is an expressly disposable product and cannot be reused. Dial medicali srl is responsible only for the safety and efficacy in case of first use. Any reuse of the device may cause contamination risks. For disposal, comply with the regulations in force.

CONTRAINDICATIONS

Contraindications are the following, but can be evaluated by the surgeon:

- Presence of infection
- Patients with mental instability, such as not to allow to follow the instructions provided by the surgeon for the post-operative course
- Patient HIVpositive
- Patients with severe osteoporosis
- Patient with diabetes mellitus
- Patients who have shown a particular sensitivity or allergy to the implant material of the PinFix® systems

IMPORTANT NOTES

- PinFix® has two configurations: cylindrical and conical
- When you use cylindrical pin you can retract it from the bone without losing its seal
- When you use conical pin never retract it after its insertion to avoid mechanical loosening
- Otherwise the conical pin can be screwed deeper to increase grip in the bone during treatment.
- Not all surgeries are successful. Further complications can develop at any time due to misuse, medical reasons, or device failures, resulting in the need for new surgery to replace or remove the device. Pre-operative and operative procedures, which require knowledge of surgical techniques, the correct choice and positioning of the devices are fundamental factors for the success of the surgery. Correct information to the patient on the post-operative course is an important factor for the success of the surgery.

WARNINGS

- In pediatric applications pay attention to the positioning of the screw in the metaepiphyseal area to prevent possible interference with the growth cartilages (physis).
- Select the diameter and length of the screw based on the size of the bone; the PinFix® screw with D. 6 mm should be used on bone with a diameter of at least 18/20 mm, otherwise evaluated to use a smaller caliber screws.
- If the screw is inserted by hand, ushe dedicated key or spindle-Jacobs type.
- If the screw is inserted by a motor, do not use a high speed rotation.
- The insertion of the pin-wire must take place by rotation of the same using hand or electric power and not by percussio
- To assure the maximum hold the screw has to always bebi-cortex with at least 1-2 threads showing. Not respecting this procedure reduces the mechanical stability of the screw (around 30%-50%) Moreover, it can cause premature osteolysis.
- To prevent osteolysis and to even out the load between the bone and the fixator the screws should be positioned perpendicular to the diaphysis or parallel the articular joint.
- Self-drilling screws do not normally require pre-drilling but the same may be necessary if the cortical bone is particularly thick; if necessary, use 4.5 / 4.8mm perforators for 6mm diameter screws and 3.2 / 3.5mm perforators for 5mm diameter screws; if necessary, screws with a diameter of 3 - 4 mm require a perforator with a diameter of 1.0 / 1.2mm smaller than the gauge of the screw used.
- Stainless steel PinFix® screw contains Nickel. It is recommended that the patient undergoes nickel allergy tests if there is any doubt.
- The PinFix® device has not been evaluated for safety and compatibility in the magnetic resonance (MR) environment and has not been tested for heating or migration to the MR environment.
- It is mandatory to inform patients about the hygiene to be carried out at the screw application site
- When removing the screws, especially if coated with hydroxyapatite, the surgeon is required to evaluate any patient sedation to prevent any discomfort / pain sensation in the patient himself
- Pay particular attention not to penetrate the joint with the screws

ADVERSE REACTIONS

- Damage to nerves and blood vessels
- Bone infection
- Edema, swelling and compartment syndrome
- Dislocation, subluxation, joint stiffness
- Fractures caused by the screw holes themselves, following their removal
- Breakage of the screw
- Tissue necrosis
- Post-operative bleeding or vascular disorders

MATERIALS, BRANDS AND PRODUCTION

PinFix® device is manufactured according to ISO 5832-1 / 3 standards (316L stainless steel or Ti6Al4V titanium), PinFix®HA has thread with hydroxyapatite (Osprovit® by LincoTek SpA).

PinFix® is produced by DialMedicali srl Milano (Italy) and also marketed with its registered trademark Dialmed®

CE0476

The packaged contains a closing cap.

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