




*Follow  
the future...*





Trauma struts for an easy reduction

CARBON FIBER RINGS or ALLUMINIUM RINGS

A new fixator which permits a flexible or symmetrical mounting

CHOOSE THE BEST SOLUTION




No spheres are needed any more

**Fixed position**

For a faster pre-mounting, the white circles will help you


Align the white "P" circle with the "D" hole in the distal ring

Put the same length to the same struts  
Remember: you can skip to flexible mounting whenever you like. Just move the strut to the desired hole.



**Extension plate**

You can use an extension plate for rings very close to each other  
You can choose to add extension plate either during the theatre or during the treatment.  
You can add this information in the software quite easily.  
P.S.: Look in the software updates to know where to put this information.



## 1 RED or BLACK rings

After filling the patient data, the software will ask you to select the color of the rings: RED (I-Fix<sup>1</sup>) or BLACK (I-Fix<sup>2</sup>)  
Remember you can't mix the colour of the rings during the surgery.

**Treatment**

Angulations:		Translations:	
Valgus: 0° over 18.91°		Medial: 0mm over 5.64mm	
Procurvatum: 0° over 2.84°		Anterior: 0mm over 3.5mm	
External: 0° over 0°		Short: 4mm over 5mm	

Date	Weekday	211 L	140 M	209 L	152 M	226 XL	211 L	TOTAL	Xtension Plate
25/05/23	Thu							TOTAL	
26/05/23	Fri	211 L	140 M	209 L	152 M	226 XL	211 L	TOTAL	
27/05/23	Sat	212 L	141 M	209 L/XL	153 M	226 XL	212 L	TOTAL	
28/05/23	Sun	213 L	142 M	210 L	154 M/L	227 XL	213 L	TOTAL	
29/05/23	Mon	214 L	143 M	210 L	155 M/L	228 XL	214 L	TOTAL	

End of axial translation



You can discover which are the corrections to be done for each day



## 2 XTENSION PLATE

You can put the xtension plate either during OR or during the treatment.

If you decide to put the xtension plate in OR, the software at page n°2 "Mounting Parameters" you will see the image of the xtension plate on each strut, click on the button in order to choose the place you mount it.



You can put the xtension plate either during OR or during the treatment.

If you decide to put the xtension plate during the treatment, go to prescription page, click on the xtension plate image in order to choose the correct position of the xtension plate.

